

Small Variants Summary

	SNVs	Insertions	Deletions
Total Passing	3,564,342	280,773	292,298
Percent Found in dbSNP	99.18%	80.29%	82.87%
Het/Hom Ratio	1.62	1.91	2.21
Ts/Tv Ratio	2.08	-	-

Structural Variants Summary

SV Type	Total	Number in Genes
CNV	89	54
Insertion	686	343
Tandem Duplication	139	84
Deletion	4,309	2,290
Inversion	127	61

Coverage Histogram




Figure 3: Resequencing Report— Variant tables and a coverage histogram from the HiSeq Analysis Software v2.0 WGS Analysis Pipeline are displayed to assess the quality of each genome analyzed.

Materials and Methods

Library Preparation and Sequencing

Sequencing libraries for human WGS were prepared from sample NA12878 (Coriell Institute for Medical Research) using the TruSeq® DNA PCR Free Library Preparation Kit (Illumina, Catalog No. FC-121-3001). The libraries were sequenced to 40x on a HiSeq System, using paired end 2 × 100 read length.

For paired tumor-normal analysis, sequencing libraries were prepared from the cell line HCC2218 (University of Texas Southwestern Medical Center). For large somatic variants, paired tumor-normal samples were prepared from HCC1187 (University of Texas Southwestern Medical Center). The libraries were prepared using the TruSeq DNA PCR Free Library Preparation Kit and sequenced on a HiSeq System, using paired end 2 × 100 read length. The paired tumor-normal samples were sequenced to 80x and 40x respectively.

For Research Use Only. Not for use in diagnostic procedures.

Somatic Small Variants Summary

	SNVs	Deletions	Insertions
Total	15,345	551	421
Number in Genes	6,785	246	200
Number in Exons	284	13	10
Number in Coding Regions	169	8	4
Splice Site Region	16	0	1
Stop Gained	8	0	0
Stop Lost	0	0	0
Frameshift	0	6	4
Non-synonymous	130	2	0
Synonymous	31	0	0
Mature miRNA	0	0	0
UTR Region	115	5	6
dbSNP	1,608	16	41

Somatic Structural Variants Summary

SV Type	Total	Number in Genes
CNV	169	162
Deletion	156	109
Tandem Duplication	93	85
Insertion	0	0
Inversion	19	16
Translocation Breakend	38	19

Depth/B-Allele Plot




Figure 4: Somatic Analysis Report— Tables and graphics from HiSeq Analysis Software v2.0 Tumor Normal Analysis Pipeline summarizing somatic variant analysis of a tumor-normal sample pair.

Whole-Genome Sequencing Pipelines

WGS data were analyzed with the HiSeq Analysis Software v2.0 WGS Analysis Pipeline, which includes alignment with the Isaac Aligner and identification of SNVs and small indel calls with Starling Small Variant Caller. For large variant calling, the Canvas CNV Caller predicted copy number deletions and duplications, while deletions, insertions, and tandem duplications were predicted with the Manta Structural Variant Caller. For quality comparison, WGS data was also analyzed with the BWA/GATK v3.0 pipeline.

