

Fragmentation Master Mix.

Illumina Stranded mRNA Prep, Ligation Checklist

For Research Use Only. Not for use in diagnostic procedures.

| Purify and Fragment mRNA | Nuclease-free ultrapure water (10.5 μl) | 42 Place on the magnetic stand and wait 2 minutes. |
|--|--|--|
| ☐ 1 Dilute 25–1000 ng total RNA in nuclease-free | ► EPH3 (10.5 μl)□ 21 Centrifuge the sealed PCR plate at 280 × g for | □ 43 Transfer 17 µl supernatant. |
| ultrapure water to 25 µl. | 10 seconds. | 44 Set aside on ice. |
| 2 Add 25 µl RPB. | □ 22 Add 25 µl BBB to each well. | E 44 Oct dolde office. |
| 3 Mix using either method: | 23 Mix using either method: | |
| Shake at 2000 rpm for 1 minute, and then | Shake at 2000 rpm for 1 minute, and then | |
| centrifuge at 280 × g for 10 seconds. | centrifuge at 280 × g for 10 seconds. | |
| ▶ Pipette 10 times. | Pipette 10 times. | |
| 4 Place on the thermal cycler and run the | 24 Incubate at room temperature for 5 minutes. | |
| mRNA_CAP program. | 25 Place on the magnetic stand and wait 2 | |
| \square 5 Centrifuge at 280 × g for 10 seconds. | minutes. | |
| \square 6 Place on the magnetic stand and wait 2 | \square 26 Remove and discard 50 μ l supernatant. | |
| minutes. | 27 Remove from the magnetic stand. | |
| 7 Remove and discard supernatant. | □ 28 Add 100 µl BWB. | |
| 8 Remove from the magnetic stand. | 29 Mix using either method: | |
| □9 Add 100 µl BWB. | Shake at 2000 rpm for 1 minute, and then | |
| ☐ 10 Mix using either method: | centrifuge at $280 \times g$ for 10 seconds. | |
| Shake at 2000 rpm for 1 minute, and then | Pipette 10 times. | |
| centrifuge at $280 \times g$ for 10 seconds. | \square 30 Place on the magnetic stand and wait 2 | |
| Pipette 10 times. | minutes. | |
| ☐ 11 Place on the magnetic stand and wait 2 | \square 31 Remove and discard supernatant. | |
| minutes. | ☐ 32 Remove residual BWB. | |
| 12 Remove and discard supernatant. | \square 33 Remove from the magnetic stand. | |
| 13 With a 20 µl pipette, remove all residual BWB. | 34 Thoroughly pipette Fragmentation Master Mix. | |
| 14 Remove from the magnetic stand. | ☐ 35 Add 19 µl Fragmentation Master Mix. | |
| 15 Add 25 µl ELB. | ☐ 36 Mix using either method: | |
| 16 Mix using either method: | Shake at 2200 rpm for 1 minute. | |
| Shake at 2200 rpm for 1 minute. | Pipette until resuspended. | |
| Pipette until resuspended. | \square 37 If shaking did not fully resuspend, pipette until | |
| 17 If shaking did not fully resuspend, pipette until | resuspended. | |
| resuspended. | 38 Incubate at room temperature for 2 minutes. | |
| 18 Centrifuge at 280 × g for 10 seconds. | ☐ 39 Centrifuge at 280 × g for 10 seconds. | |
| 19 Place on the thermal cycler and run the | 40 Place on the thermal cycler and run the | |
| mRNA_ELT program. | DEN94_8 program. | |
| ☐ 20 Combine the following volumes to prepare | \square 41 Centrifuge at 280 \times g for 10 seconds. | |

Illumina Stranded mRNA Prep, Ligation Checklist

Synthesize First Strand cDNA

| □ 1 | Combine the following volumes to prepare |
|-------------|--|
| | First Strand Synthesis Master Mix. |
| | FSA (9 μl) |
| | P RVT (1 μl) |
| \square 2 | Thoroughly pipette First Strand Synthesis |
| | Master Mix. |
| \square 3 | Centrifuge at $280 \times g$ for 10 seconds. |
| $\Box 4$ | Add 8 µl First Strand Synthesis Master Mix. |
| □ 5 | Pipette 10 times. |
| □ 6 | Place on the thermal cycler and run the FSS |
| | program. |
| | - |

Synthesize Second Strand cDNA

| □ 1 | Centrifuge at 280 × g for 10 seconds. |
|--------------|---|
| \square 2 | Add 25 µl SMM. |
| \square 3 | Pipette 10 times. |
| \Box 4 | Place on the thermal cycler and run the |
| | SSS program. |
| \Box 5 | Centrifuge at 280 × g for 10 seconds. |
| □6 | Add 90 µl AMPure XP. |
| \Box 7 | Mix using either method: |
|) | Shake at 2000 rpm for 1 minute, and then |
| | centrifuge at $280 \times g$ for 10 seconds. |
|) | Pipette until resuspended. |
| 8 | Incubate at room temperature for 5 minutes. |
| 9 | Place on the magnetic stand and wait 5 |
| | minutes. |
| | Remove and discard 130 µl supernatant. |
| | Wash beads as follows. |
| | a Add 175 µl fresh 80% EtOH. |
| | b Wait 30 seconds. |
| | c Remove and discard supernatant. |
| | Repeat wash a second time. |
| | Remove residual EtOH. |
| | Air-dry for 2 minutes. |
| | Remove from the magnetic stand. |
| | Add 19.5 µl RSB. |
| □1/ | Mix using either method: |
|) | Shake at 2200 rpm for 1 minute. |
| □ 40 | Pipette until resuspended. |
| □18 | If shaking did not fully resuspend, pipette until |
| | resuspended. |
| | Incubate at room temperature for 2 minutes. |
| | Centrifuge at 280 × g for 10 seconds. |
| □ 21 | Place on the magnetic stand and wait 2 |
| | minutes. |
| \square 22 | Transfer 17.5 µl supernatant. |

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

Illumina Stranded mRNA Prep, Ligation Checklist

For Research Use Only. Not for use in diagnostic procedures.

Adenylate 3' Ends

| □ 1 Add 12.5 µl ATI | _4 |
|---------------------|----|
|---------------------|----|

- □2 Using a 200 µl pipette, pipette 10 times.
- Place on the thermal cycler and run the ATAIL program.

Ligate Anchors

- \square 1 Centrifuge at 280 × g for 10 seconds.
- 2 Add the following volumes *in the order listed*.

| Reagent | Volume for Input ≤ 100 ng (µI) | Volume for Input > 100 ng (μΙ) |
|----------------------|--------------------------------------|--------------------------------------|
| RSB | 2.5 | 0 |
| RNA Index Anchors | 2.5 | 5 |
| LIGX | 2.5 | 2.5 |

- ☐3 Pipette 10 times, and then seal.
- Place on the thermal cycler and run the LIG program.
- \square 5 Centrifuge at 280 × g for 10 seconds.
- \square 6 Add 5 μ I STL.
- \square 7 Pipette 15 times to mix.

Clean Up Fragments

- ☐ 1 Add 34 µl AMPure XP.
- ☐ 2 Mix using either method:
 - Shake at 2000 rpm for 1 minute, and then centrifuge at 280 × g for 10 seconds.
 - Pipette until resuspended.
- □3 Incubate at room temperature for 5 minutes.
- ☐ 4 Place on the magnetic stand and wait 5 minutes.
- □ 5 Remove and discard 67 μl supernatant.
- ☐ 6 Wash beads as follows.
 - □a Add 175 µl fresh 80% EtOH.
 - ☐ b Wait 30 seconds.
 - C Remove and discard supernatant.
- \square 7 Wash beads a **second** time.
- 8 Remove residual EtOH.
- \square 9 Air-dry for 2 minutes.
- □ 10 Remove from the magnetic stand.
- □ 11 Add 22 µl RSB.
- ☐ 12 Mix using either method:
 - Shake at 2200 rpm for 1 minute.
 - Pipette until resuspended.
- 13 If shaking did not resuspend, pipette until resuspended.
- ☐ 14 Incubate at room temperature for 2 minutes.
- \square 15 Centrifuge at 280 × g for 10 seconds.
- ☐ 16 Place on the magnetic stand and wait 2 minutes.
- □ 17 Transfer 20 µl supernatant

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.



Illumina Stranded mRNA Prep, Ligation Checklist

Clean Up Library

For Research Use Only. Not for use in diagnostic procedures.

| Amplify Library |
|--|
| □ 1 Pierce the index adapter plate wells. □ 2 Add the following volumes in the order listed. ▶ UDP0XXX (10 µl) ▶ EPM (20 µl) □ 3 Pipette 10 times. □ 4 Place on the thermal cycler and run the PCR program. |
| SAFE STOPPING POINT |
| If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days. |

| □ 1 Centrifuge at 280 × g for 10 seconds. □ 2 Add 50 µl AMPure XP. |
|--|
| 3 Mix using either method: |
| Shake at 2000 rpm for 1 minute, and ther centrifuge at 280 × g for 10 seconds. Pipette until resuspended. |
| 4 Incubate at room temperature for 5 minutes. |
| 5 Place on the magnetic stand and wait 5 minutes. |
| 6 Remove and discard 90 µl supernatant. |
| 7 Wash beads as follows. |
| a Add 175 µl fresh 80% EtOH. |
| □b Wait 30 seconds. |
| \square c Remove and discard supernatant. |
| 8 Wash beads a second time. |
| 9 Remove residual EtOH. |
| □ 10 Air-dry for 2 minutes. |
| \square 11 Remove from the magnetic stand. |
| □ 12 Add 17 µl RSB to each well. |
| ☐ 13 Mix using either method: |
| Shake at 2200 rpm for 1 minute. |
| Pipette until resuspended. |
| □ 14 If shaking did not resuspend, pipette until |
| resuspended. |
| \square 15 Incubate at room temperature for 2 minutes. |
| \square 16 Centrifuge at 280 \times g for 10 seconds. |
| ☐ 17 Place on the magnetic stand and wait 2 |
| minutes. |
| □ 18 Transfer 15 µl supernatant. |
| SAFE STOPPING POINT |

Check Library

Analyze 1 µl library using the Agilent 2100
 Bioanalyzer and DNA 1000 Kit.
 [Optional] Analyze 2 µl library using the Qubit dsDNA BR Assay Kit.

If you are stopping, seal the plate and store at

-25°C to -15°C for up to 7 days.

Dilute Library to the Starting Concentration

- ☐ 1 Obtain the molarity value:
 - ▶ Bioanalyzer quantification only—Use the molarity value obtained for the library.
 - Bioanalyzer and Qubit quantification— Calculate molarity value using the average size and concentration.
- Using the molarity value, calculate the volumes of RSB and library needed to dilute libraries to the starting concentration for your system.

| Sequencing System | Starting Concentration (nM) | Final Loading Concentration (pM) |
|-----------------------------------|-----------------------------------|--|
| NextSeq 550 and NextSeq 500 | 1 | 1.1–1.4 |
| NovaSeq 6000 | 0.5 | 100 |

- Dilute each library to the starting concentration. Combine 10 µl each diluted library in a tube.
- 4 Follow denature and dilute instructions to dilute libraries.

Acronyms

| Acronym | Definition |
|---------|---|
| ATL4 | A-Tailing Mix |
| BBB | Bead Binding Buffer |
| BWB | Bead Washing Buffer |
| cDNA | Complementary DNA |
| ELB | Elution Buffer |
| EPH3 | Elute, Prime, Fragment High Concentration Mix |
| EPM | Enhanced PCR Mix |
| EtOH | Ethanol |
| FSA | First Strand Synthesis Act DMix |
| LIGX | Ligation Mix |
| mRNA | Messenger RNA |
| RPBX | RNA Purification Beads |
| RSB | Resuspension Buffer |
| SMM | Second Strand Master Mix |
| STL | Stop Ligation Buffer |
| | |