


Illumina 5-Base DNA Prep with Enrichment

A single assay for simultaneous targeted detection of genomic variants and methylation events



High-sensitivity detection of methylation events using novel chemistry



Two-in-one assay features a streamlined workflow and easy analysis



Flexible customization for targeted genome and methylome insights

High-sensitivity multiomic detection

DNA is inherently multiomic, holding both genetic and epigenetic molecular information. Beyond the sequence of adenine (A), thymine (T), guanine (G), and cytosine (C), there are modified bases such as 5-methylcytosine (5mC) that help direct gene expression (Figure 1). Detecting both genomic variation and DNA methylation can reveal hidden mechanisms of health and disease. Studying the genome and methylome typically requires separate assays and data analysis steps. Illumina 5-Base DNA Prep with Enrichment leverages a novel enzymatic method and optimized analysis for detection of five bases (A, T, G, C, and 5mC) from the same sample, in the same assay, using the same sequencing reads. With hybrid-capture enrichment and custom-designed panels, users can focus on specific targets of interest, enabling cost-effective, deep coverage with high-sensitivity multiomic detection.

Streamlined, flexible workflow

Illumina 5-Base DNA Prep with Enrichment enables targeted DNA sequencing and methylation sequencing with one easy-to-use assay (Figure 2). This single-vendor solution provides a streamlined library-to-analysis workflow that can be completed in less than three days. Illumina 5-Base DNA Prep with Enrichment is compatible with cell-free DNA (cfDNA) and genomic DNA (gDNA) from blood, cell lines, or fresh frozen tissue. Optimized enrichment library prep, which includes a simple, one-step 5mC-to-T base conversion, requires minimal touchpoints and is completed in fewer than 11 hours (Table 1, Figure 3, Figure 4).^{*} Illumina 5-Base DNA Prep with Enrichment uses custom hybrid-capture probe panels for targeted studies. The assay supports four-plex enrichment and up to 384 unique dual indexes. Sequence libraries on the NextSeq™ 2000 System, NovaSeq™ 6000 System, NovaSeq 6000Dx Instrument (RUO mode), or the NovaSeq X Series, followed by simplified data analysis. An integrated DRAGEN™ pipeline generates a dual readout (Figure 5, Figure 6).

^{*} Based on processing 24 samples manually; includes DNA shearing (if applicable), and library quantification. Normalization and pooling times are not included.




Figure 1: DNA methylation of C to 5mC is a well-studied epigenetic mark for gene regulation

Illumina 5-Base DNA Prep with Enrichment detects 5mC along with unmodified A, T, G, and C bases, providing both genomic and epigenomic insights from a single NGS assay.





Table 1: Illumina 5-Base DNA Prep with Enrichment parameters


	Input	Library prep time ^a
Genomic DNA	50–100 ng	< 11 hr
Cell-free DNA	1–20 ng	< 11 hr

^a. Based on processing 24 samples manually; includes DNA shearing (if applicable), and library quantification. Normalization and pooling times are not included.

Novel chemistry for direct conversion of 5-methylcytosine to thymine

Illumina 5-Base DNA Prep with Enrichment is a fundamentally different approach to targeted genome and methylome analysis. Traditional methods for detecting DNA methylation use bisulfite treatment or enzymes to convert unmethylated cytosine to thymine (Figure 3B). Because most cytosines in the genome are unmodified, this approach greatly reduces nucleotide diversity, making sequencing reads harder to align. Illumina 5-Base DNA Prep with Enrichment uses a novel engineered enzyme to directly convert only 5mC to thymine in a single incubation step during library preparation (Figure 3A). Illumina 5-base chemistry preserves library complexity and ensures high data alignment rate, increasing the sensitivity of detection.





Single readout with combined genome and methylome data

Integrated DRAGEN™ secondary analysis provides accurate annotation of both methylation and genomic variants in a single readout (Figure 5, Figure 6). Novel methylation-aware algorithms leverage the complementary strand sequence to discern between a thymine indicating a methylation event and a thymine that represents a single nucleotide variant (SNV) (Figure 5). Integrated unique molecular identifiers (UMIs) help provide confident single-molecule resolution of methylation status. Duplex UMI collapsing improves technical reproducibility for detecting small methylome differences and allows sensitive interrogation of asymmetric methylation across cytosine-guanine dinucleotide (CpG) sites (Figure 5).

Illumina 5-base methylation reporting features are available within the DRAGEN Enrichment pipeline with an easy checkbox option. Secondary analysis can be performed via BaseSpace™ Sequence Hub or Illumina Connected Analytics cloud platforms, or on a DRAGEN server.





Figure 5: DRAGEN algorithms distinguish methylation events from SNVs

Innovative DRAGEN algorithms leverage the complementary strand sequence to accurately discern between methylation and small variant calls in the same read. For 5mC converted to T, the complementary base will be G, whereas for a C-to-T genomic variant, the complementary base will be A. Duplex UMI collapsing further enables single-molecule resolution to detect symmetric and asymmetric methylation at CpG sites.



Customization to focus on regions of interest

Illumina 5-Base DNA Prep with Enrichment requires panels specifically designed for the 5-base assay. Panel design is key for optimal performance and efficient capture of both methylated and nonmethylated targets.[†] To focus sequencing on regions of interest, Illumina Custom Enrichment Panel v2, in a broad range of panel sizes, can be designed and ordered using DesignStudio™ software. The specialized design algorithm accommodates the 5mC-to-T base conversion unique to the Illumina 5-base solution.[†] DesignStudio software supports design of custom panels from target gene lists or genomic coordinates. Content from existing panels can be used as a starting point. The DesignStudio tool also allows users to select regulatory regions, such as promoters, within defined targets.

High-sensitivity detection

Illumina 5-Base DNA Prep with Enrichment is optimized to meet sensitive detection requirements with exceptional performance for multiple applications. Detect low allele frequency variants from low-input samples with integrated UMIs. With complementary workflows, users can scan the whole-genome and methylome using Illumina 5-Base DNA Prep for broad discovery, then use those results to guide targeted panel design for deeper sequencing using Illumina 5-Base DNA Prep with Enrichment.

[†]Standard enrichment panels do not account for methylated and nonmethylated bases. Panels made for methods that convert unmethylated C, such as bisulfite sequencing, are designed for 3-base nucleotide diversity (Figure 3) and are not compatible with Illumina 5-Base DNA Prep with Enrichment.

Summary

Combining analysis of genetic variation and DNA methylation helps maximize insights from every sequencing read. Illumina 5-Base DNA Prep with Enrichment offers a library-to-analysis solution for simultaneous targeted genome and methylome profiling in one optimized workflow. Custom enrichment panel design focuses sequencing on regions of interest. Novel chemistry and analytical techniques deliver high sensitivity of detection of genomic variants and methylation in a single assay.

Learn more →

[Illumina 5-Base DNA Prep with Enrichment](#)



1.800.809.4566 toll-free (US) | +1.858.202.4566 tel
techsupport@illumina.com | www.illumina.com

© 2025 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.
M-GL-03691 v1.0

Ordering information

Product	Catalog no.
Library prep	
Illumina 5-Base DNA Prep with Enrichment (24 samples)	20140366
Illumina 5-Base DNA Prep with Enrichment (96 samples)	Coming soon
Panels	
Illumina Custom Enrichment Panel v2 (32 µl, 120 bp)	20073953
Illumina Custom Enrichment Panel v2 (384 µl, 120 bp)	20073952
Illumina Custom Enrichment Panel v2 (1536 µl, 120 bp)	20111339
Indexes	
Illumina DNA/RNA UD Indexes Set A, Tagmentation (96 indexes, 96 samples)	20091654
Illumina DNA/RNA UD Indexes Set B, Tagmentation (96 indexes, 96 samples)	20091656
Illumina DNA/RNA UD Indexes Set C, Tagmentation (96 indexes, 96 samples)	20091658
Illumina DNA/RNA UD Indexes Set D, Tagmentation(96 indexes, 96 samples)	20091660
Illumina Unique Dual Indexes, LT (48 indexes, 48 samples)	20098166
Analysis	
Illumina DRAGEN server v4	20051343
Illumina Analytics - 1 iCredit	20042038
Illumina Analytics Starter Package - 1000 iCredits	20042039
Illumina Analytics - 5000 iCredits	20042040
Illumina Analytics - 50,000 iCredits	20042041
Illumina Analytics - 100,000 iCredits	20042042