

Rapid Viral Surveillance Panel v2

Streamlined whole-genome sequencing for high-risk viral surveillance and research



Expanded panel provides coverage of ~200 viruses, including those of public health concern¹⁻⁶



Hybrid-capture enrichment accommodates RNA and DNA viral pathogens



Rapid workflow supports a range of host and environmental sample types⁷

Identifying high-impact viruses for public health surveillance

Viral genome surveillance plays a pivotal role in global health security by providing invaluable insights into the evolution, spread, and activity of pathogens.¹ Analyzing the genetic makeup of viruses using next-generation sequencing (NGS) allows scientists to track mutations that may affect transmissibility, virulence, or resistance to treatment.⁸ This information is crucial for designing effective diagnostic tests, therapeutics, and vaccines to combat emerging infectious diseases.

The Illumina Rapid Viral Surveillance Panel v2 is an NGS panel that enables the detection and whole-genome sequencing (WGS) of ~200 viruses, including viruses identified as important risks to public health (Table 1).¹⁻⁶ The panel uses a hybrid-capture target enrichment workflow that allows for sequencing of various sample types in as little as 24 hours without the high read depth required by shotgun metagenomics methods. Compared to other targeted resequencing methods, such as

amplicon sequencing, hybrid-capture provides more uniform coverage across viral genomes and a greater ability to identify mutations and divergent sequences, making the Viral Surveillance Panel v2 ideal for outbreak surveillance and variant monitoring.

Streamlined NGS workflow

The Rapid Viral Surveillance Panel v2 workflow enriches for viral genomes from a range of sample types, including wastewater, serum, plasma, skin lesions, and nasopharyngeal swabs.⁷ Libraries are prepared from RNA or DNA extracted from host or environmental samples, sequenced on an Illumina benchtop sequencing system, and analyzed using the DRAGEN™ Microbial Enrichment Plus App available on BaseSpace™ Sequence Hub. The library preparation and sequencing steps can be completed in as little as 24 hours with minimal hands-on time (Figure 1).⁷




Figure 1: Rapid Viral Surveillance Panel v2 workflow

In a rapid, streamlined, and comprehensive workflow, libraries are prepared from environmental or host samples, sequenced on an Illumina sequencing system, and analyzed with the DRAGEN Microbial Enrichment Plus App for viral detection, whole-genome consensus generation, read mapping to viral best hits, and strain typing. Sequencing can be completed in under two days, but the exact time varies with sample read depth and sequencing system used.

- a. Sequencing time for 25M-100M fragments at 2 × 101 bp is ~7-8 hr on the MiSeq i100 Series.
- b. Sample prep and sequencing steps can be completed in less than one day using the MiSeq i100 Series. Sequencing time varies based on the sequencing system and read depth.

Table 1: Key high-risk viruses included on the Rapid Viral Surveillance Panel v2

| | | | |
|---------------------------------------|------------------------------------|--|---|
| Adeno-associated virus 2 | Human adenovirus A–G | Mayaro virus | Sabia virus |
| Aichi virus 1 | Human bocavirus | Measles virus | Salivirus A |
| Aigai virus | Human coronavirus | Menangle virus | Sandfly fever Sicilian virus |
| Bombali virus | Human cytomegalovirus | Middle East respiratory syndrome–related coronavirus | Sapovirus |
| Bourbon virus | Human immunodeficiency virus 1/2 | Mpox virus | Severe acute respiratory syndrome–related coronavirus |
| Cache Valley virus | Human metapneumovirus | Mumps virus | Severe acute respiratory syndrome–related coronavirus 2 |
| California encephalitis virus | Human papillomavirus | Murray Valley encephalitis virus | Semliki Forest virus |
| Chapare virus | Human parainfluenza virus 1–4 | Nipah virus | Severe fever with thrombocytopenia syndrome virus |
| Chikungunya virus | Human parechovirus | Norovirus | Sindbis virus |
| Colorado tick fever virus | Human parvovirus B19 | Omsk hemorrhagic fever virus | Snowshoe hare virus |
| Coxsackievirus A/B | Influenza virus A–C | O'nyong-nyong virus | Sosuga virus |
| Crimean-Congo hemorrhagic fever virus | Jamestown Canyon virus | Oropouche virus | St. Louis encephalitis virus |
| Dengue virus 1–4 | Japanese encephalitis virus | Poliovirus | Tacheng tick virus 2 |
| Ebola virus | Junin virus | Polyomavirus | Tahyna virus |
| Echovirus | Kyasanur Forest disease virus | Powassan virus | Tick-borne encephalitis virus |
| Enterovirus A–D | La Crosse virus | Punta Toro virus | Torque teno virus |
| Epstein-Barr virus | Lassa virus | Rabies virus | Toscana virus |
| Equine encephalitis virus | Lloviu virus | Ravn virus | Usutu virus |
| Guanarito virus | Lujo virus | Respiratory syncytial virus A/B | Varicella-zoster virus |
| Hantavirus | Lymphocytic choriomeningitis virus | Rhinovirus A–C | Variola virus |
| Heartland virus | Lyssavirus | Rift Valley fever virus | West Nile virus |
| Henipavirus | Machupo virus | Ross River virus | Yellow fever virus |
| Hepatovirus A–E | Mamastrovirus | Rotavirus A/B/C/H | Zika virus |
| Herpes simplex virus 1/2 | Marburg virus | Rubella virus | |

Library preparation

The Rapid Viral Surveillance Panel v2 library preparation workflow consists of pre-enrichment and enrichment steps. The rapid pre-enrichment library step is greatly simplified and streamlined to be performed in half the time, with a reduced number of steps, compared to the nonrapid format. Pre-enrichment generates hundreds of thousands of nontargeted libraries that are enriched with Viral Surveillance Panel v2 probes using a hybrid-capture approach. Enrichment with on-bead tagmentation provides a rapid, automation-compatible workflow that can be completed with minimal hands-on time. The protocol accommodates sample input amounts ranging from 10 ng to 100 ng total nucleic acid and supports multiplexing of up to 384 samples in a single run.

Sequencing

The lower read depth requirements for libraries enriched with Viral Surveillance Panel v2 allow for multiple sequencing system options, including the benchtop MiniSeq™, MiSeq™ i100 Series, NextSeq™ 550, NextSeq 1000, and NextSeq 2000 Systems. Viral titer, nucleic acid sample quality, sample read depth, and the number of reads per sample impact the number of virus-specific reads and sequence coverage obtained. The general sequencing read depth recommendation for good quality samples is a minimum of 1M fragments per sample with a read length of 2 × 150 bp. The recommended sample read depth also varies with sample type. For more complex samples, such as wastewater, a minimum of 4M fragments per sample are recommended. Abundant off-target reads are expected if other microbial nucleic acids are present, such as from bacteria found in complex sample types.

Data analysis

Data generated using the Viral Surveillance Panel v2 are analyzed using the DRAGEN Microbial Enrichment Plus App available on BaseSpace Sequence Hub and onboard MiSeq i100 Series Systems. This easy-to-use analysis pipeline provides sample quality control, reference-guided alignment to a broad, curated viral genome database, variant calling, viral genome consensus sequence generation, antiviral resistance prediction for Influenza A/B viruses, flexible reporting options, and integration with Pangolin and Nextclade for further phylogenetic assignment of supported viruses.

Performance

Target enrichment

Compared to shotgun metagenomic sequencing, where all RNA or DNA is sequenced, targeted hybrid-capture using the Viral Surveillance Panel v2 minimizes unnecessary sequencing of host and nontargeted microbes, reducing costs and allowing for broad sequencing of viral genomes on benchtop sequencing systems.⁷

Viral genome recovery using enrichment with the Viral Surveillance Panel v2 was compared with shotgun metagenomic sequencing without enrichment. Analysis was performed using wastewater samples collected through collaborations with Colorado State University (n = 29) and Wisconsin State Lab of Hygiene (n = 32). The Viral Surveillance Panel v2 demonstrated superior viral genome recovery from multitarget contrived samples compared to shotgun metagenomic sequencing (Figure 2).




Figure 2: Rapid Viral Surveillance Panel v2 performance compared to shotgun metagenomic sequencing

Rapid Viral Surveillance Panel v2 (Rapid) produces higher genome coverage measured as reads per kilobase of transcript per million reads mapped (RPKM) compared to shotgun metagenomic RNA sequencing with another commercially available method. Increased coverage of targeted genomes Rapid Viral Surveillance Panel v2 ensures more full genomes are captured and delivers excellent performance for detection and surveillance. Analysis was performed using wastewater samples collected through collaborations with Colorado State University (n = 29) and Wisconsin State Lab of Hygiene (n = 32).

Wastewater surveillance

Surveillance for viral sequences in wastewater provides a regional indicator of communal spread of viral pathogens, giving public health professionals valuable information for response planning.⁸ The Viral Surveillance Panel can be used with these samples to enable early detection and identification of viral genomes in wastewater at lower concentrations than shotgun sequencing (Table 2).

Wastewater samples from two collection sites extracted through collaborations with the Wisconsin State Lab of Hygiene (WSLH) and Colorado State University (CSU). Three samples from each collection site were assessed. Libraries prepared from these six wastewater samples

were sequenced and normalized to 4M fragments for Viral Surveillance Panel v2 enrichment or 4M and 25M fragments for shotgun metagenomic sequencing. A sequencing depth of 4M fragments was used in this comparison because wastewater samples can vary greatly in complexity and may contain dozens of viruses at low abundance. The Viral Surveillance Panel v2 demonstrated increased viral detection sensitivity in a complex environmental sample type with low overall viral load compared to shotgun metagenomic sequencing, even when total reads for shotgun metagenomic sequencing were increased ~six-fold (Table 2).

Table 2: Top viruses detected in wastewater using Rapid Viral Surveillance Panel v2 or shotgun metagenomic sequencing

| Virus (strain) | 4M fragments | | | | 25M fragments | |
|--|-----------------------------------|------------|----------------------|------------|----------------------|------------|
| | Rapid Viral Surveillance Panel v2 | | Shotgun metagenomics | | Shotgun metagenomics | |
| | Genome coverage (%) | Read count | Genome coverage (%) | Read count | Genome coverage (%) | Read count |
| Sapovirus (GII.1) | 99.7 | 219,539 | 50.0 | 57 | 84.2 | 360 |
| Human adenovirus F (Human adenovirus 41) | 100 | 104,693 | 6.4 | 18 | 23.6 | 72 |
| Human coronavirus OC43 (HCoV_OC43) | 98.1 | 23,857 | 0 | 0 | 10.0 | 26 |
| Sapovirus (GV) | 99.6 | 10,750 | 0 | 0 | 25.3 | 19 |
| Human adenovirus E | 88.4 | 6733 | 0 | 0 | 0 | 0 |
| JC polyomavirus (JCPyV) | 99.3 | 5834 | 0 | 0 | 0 | 0 |
| Mamastrovirus 9 (MAstV9) | 99.2 | 4959 | 0 | 0 | 0 | 0 |
| Mamastrovirus 1 (MAstV1) | 98.6 | 3972 | 7.5 | 5 | 22.0 | 13 |
| Human adenovirus A (Human adenovirus 31) | 81.1 | 3449 | 0 | 0 | 0 | 0 |
| Mamastrovirus 6 (MAstV6) [MLB1] | 97.2 | 3181 | 0 | 0 | 17.3 | 9 |
| Norovirus (G1) | 96.9 | 1972 | 0 | 0 | 0 | 0 |
| BK polyomavirus (BKPyV) | 100 | 1522 | 0 | 0 | 11.1 | 4 |
| Mamastrovirus 8 (MAstV8) [VA2] | 92.1 | 1208 | 0 | 0 | 6.1 | 4 |
| Human papillomavirus 59 (HPV59; high-risk) | 69.3 | 1015 | 0 | 0 | 0 | 0 |
| Enterovirus A (not Coxsackievirus) [Enterovirus A71] | 70.0 | 295 | 5.1 | 4 | 9.0 | 6 |

Summary

The Rapid Viral Surveillance Panel v2 is part of an optimized, comprehensive workflow for detecting viral outbreaks, zoonotic surveillance, and tracking mutations. This kit uses updated chemistry to deliver the same high-quality results as the existing Viral Surveillance Panel v2 with a faster turnaround time as little as 24 hours for time sensitive applications.⁷ The streamlined workflow is compatible with a range of sample types and applications, including biological research samples and wastewater surveillance for regional presence of viruses.

The kit includes hybrid-capture probes for identifying ~200 RNA and DNA virus genomes that have been designated as high risks to public health. The hybrid-capture target enrichment minimizes the need for high sample read depth by focusing on target sequences, reducing costs while increasing throughput. Data generated using the Rapid Viral Surveillance Panel v2 can be analyzed with the user-friendly DRAGEN Microbial Enrichment Plus App onboard the MiSeq i100 Series or online through BaseSpace Sequence Hub. This robust NGS workflow delivers excellent viral capture performance for identifying DNA and RNA in complex samples, providing public health organizations and researchers with an advanced alternative to shotgun sequencing.

Learn more →

[Rapid Viral Surveillance Panel v2](#) (link pending)

[DRAGEN Microbial Enrichment Analysis Plus App](#)

[Illumina sequencing systems](#)

Ordering information

| Product | Catalog no. |
|---|-------------|
| Illumina Rapid Viral Surveillance Panel v2 Kit (96 samples) | 20158826 |

References

- Ling-Hu T, Rios-Guzman E, Lorenzo-Redondo R, Ozer EA, Hultquist JF. [Challenges and opportunities for global genomic surveillance strategies in the COVID-19 era](#). *Viruses*. 2022;14(11):2532. doi:10.3390/v14112532
- World Health Organization. Pathogens prioritization: a scientific framework for epidemic and pandemic research preparedness. [who.int/publications/m/item/pathogens-prioritization-a-scientific-framework-for-epidemic-and-pandemic-research-preparedness](#). Published July 30, 2024. Accessed March 25, 2026.
- World Health Organization. Prioritizing diseases for research and development in emergency contexts. [who.int/activities/prioritizing-diseases-for-research-and-development-in-emergency-contexts](#). Accessed March 25, 2026.
- Bloom DE, Cadarette D. [Infectious disease threats in the twenty-first century: strengthening the global response](#). *Front Immunol*. 2019;10:549. doi:10.3389/fimmu.2019.00549.
- World Health Organization. Disease Outbreak News. [who.int/emergencies/disease-outbreak-news](#). Updated March 13 2026. Accessed March 25, 2026.
- Africa Centers for Disease Control and Prevention. Diseases information. [africacdc.org/disease/](#). Accessed March 25, 2026.
- Data on file. Illumina, Inc. 2026.
- Diamond MB, Keshaviah A, Bento AI, et al. [Wastewater surveillance of pathogens can inform public health responses](#). *Nat Med*. 2022;28(10):1992-1995. doi:10.1038/s41591-022-01940-x.



1.800.809.4566 toll-free (US) | +1.858.202.4566 tel
techsupport@illumina.com | www.illumina.com

© 2026 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see [www.illumina.com/company/legal.html](#).
M-GL-04127 v1.0