

# **Correlation Engine**

## **Customer Release Notes**

**Version 2.3.0**

**February 7, 2024**

## Introduction

These Release Notes detail the latest release of Correlation Engine, including known issues.

Correlation Engine provides an interface with tools and content for exploratory, preclinical, and early translational research in drug discovery and biotechnology research. The software is a multidimensional integration platform designed to search and correlate global collections of heterogeneous, large-scale experimental data across diverse platforms, organisms, and therapeutic areas.

### NEW FEATURES:

- Images of QC plot files associated with public studies now streamed to the UI
- Updated "What's New?" section on home page
- Update copyright year to 2024

### FIXED ISSUES:

- Fixed an issue where some GWAS studies had outdated URLs in their description
- Fixed an issue where copyright was out of date

### KNOWN ISSUES:

- When searching by some biosets in Disease Atlas, some concepts show incorrect counts on the 'Studies for Concept' tab
- When searching by gene query on the cell lines tab in Body Atlas, some genes return no results in the expanded third level details

### DATA CONTENT UPDATE:

	Added since 28-Aug-23 release	Total public content
Studies	145	25,907
Biosets	363	187,883