

iScan Control Software (iCS) v4.0.5

Release Notes

For iScan Systems

April 27, 2021

FOR RESEARCH USE ONLY

© 2019 Illumina, Inc. All rights reserved.

Illumina, Infinium, iScan, , the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

Introduction

These release notes describe the key features and changes pertaining to iScan Control Software (iCS) version 4.0.5 since the last major release, version 4.0.0. This version has been tested with Windows 10 and supports both Illumina LIMS and non-Illumina LIMS customers.

I. iScan Control Software v4.0.5

NEW FEATURES:

- Enabled compatibility for iScan to upload data to Illumina Connected Analytics (ICA).

DEFECT REPAIRS:

- None

KNOWN ISSUES:

- In rare cases, Autoloader communication issues have been observed when using Illumina LIMS. This issue can be resolved by closing iCS and power cycling the instrument.

FOR RESEARCH USE ONLY

© 2019 Illumina, Inc. All rights reserved.

Illumina, Infinium, iScan, , the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.